

***Inclusive Technology***  
[www.inclusive.co.uk](http://www.inclusive.co.uk)

## **Compatibility**

The it-Control Lite is compatible with any Simply Works transmitter.

These include:

- it-Switch - a wireless switch.
- it-Send and it-Send Pro - to enable connection of wired switches.
- it-Stick - can be used as a 'wobble switch'.

## **Batteries**

The it-Control Lite requires 2 x AAA batteries. These are fitted by removing the battery cover on the underside of the device.

Ensure the batteries are inserted correctly - match the + and – symbols on the battery with those in the battery compartment. Replace the battery cover and screw. It is recommended that good quality batteries are used (such as Energiser or Duracell).

If you are not intending to use the it-Control Lite for an extended period (e.g. a number of weeks) it is advisable to remove the batteries to prevent leakage, which may damage the device.

## **Pairing it-Control Lite with a transmitter (e.g. it-Switch, it-Send)**

1. Place the it-Control Lite close to the transmitter device.
2. Press and release the PAIR button on the front of the it-Control Lite. The red light next to it will begin flashing.
3. Within 10 seconds press and release the PAIR button on the transmitter.
4. The light on the it-Control Lite will continue to flash while the devices communicate with each other.
5. Once the pairing is successful the red light on the it-Control Lite will stay on for 5 seconds.
6. If the light continues to flash for 30 seconds and then stops, repeat the process above. If it still does not pair successfully, refer to the Troubleshooting section.

Pairing information is retained in the it-Control Lite's memory even after changing the batteries.

## **Operation**

If you have paired the it-Control Lite with the it-Send or it-Switch, you must ensure the mode on these devices is set to 0.

Plug one end of the lead supplied into the socket marked OUTPUT on the it-Control Lite. Connect the other end to your preferred device. This may be an adapted battery toy or a mains control device.

The it-Control Lite activates using direct (sometimes also called 'momentary') control. This means the attached toy is activated while the switch is held down. When the switch is released the toy will stop. This is good for individuals learning about cause and effect as it provides a concrete link between the switch press and toy activation.

## **Output rating**

The output rating of the jack socket is 1 amp at up to 30V. Exceeding this rating may cause permanent damage and is not covered by the warranty.

## it-Control Lite troubleshooting

Problem	Possible cause
it-Control Lite not functioning.	<ul style="list-style-type: none"><li>• Not paired with a transmitter - see pairing instructions.</li><li>• Transmitter not set to 0 - set mode on transmitter to 0.</li><li>• No batteries in it-Control Lite, or batteries are flat – insert or replace batteries.</li></ul>
it-Control Lite works erratically or intermittently.	<ul style="list-style-type: none"><li>• Batteries in it-Control Lite flat - replace batteries.</li><li>• Out of range - move transmitter closer to it-Control Lite.</li><li>• Interference - move it-Control Lite and transmitter away from sources of interference (e.g. mobile phone).</li></ul>
Devices will not pair.	<ul style="list-style-type: none"><li>• No batteries in it-Control Lite, or batteries are flat - insert or replace batteries.</li><li>• No batteries in transmitter, or batteries are flat - insert or replace batteries.</li><li>• Ensure it-Control Lite is close to transmitter while pairing (within 1 metre).</li></ul>

## **Maintenance**

Your Simply Works device contains no user serviceable parts. If it does not work as expected consult the troubleshooting guide. If there are still problems with the device please contact your supplier.

## **Warranty**


Your Simply Works device is warranted against defects in manufacture or component failure. The unit is designed for domestic and educational applications. Use outside these areas will invalidate the warranty. Unauthorised repair or modification, mechanical abuse, immersion in any liquid or connection to incompatible equipment will also invalidate the warranty.

## **Technical support**

If you need further support with any of the Simply Works range, please contact us:

Inclusive Technology Ltd  
Riverside Court, Huddersfield Road  
Delph  
Oldham OL3 5FZ  
United Kingdom  
Tel: +44 (0)1457 819790  
Fax: +44 (0)1457 819799  
email: [inclusive@inclusive.co.uk](mailto:inclusive@inclusive.co.uk)  
[www.inclusive.co.uk](http://www.inclusive.co.uk)

## Notes



## **Inclusive Technology Ltd**

Riverside Court, Huddersfield Road  
Delph, Oldham OL3 5FZ

United Kingdom

Tel: +44 (0)1457 819790

Fax: +44 (0)1457 819799

email: [inclusive@inclusive.co.uk](mailto:inclusive@inclusive.co.uk)

[www.inclusive.co.uk](http://www.inclusive.co.uk)